

ABON™ HIV Tri-Line Rapid Tests

A Preferred Choice in HIV Rapid Tests

- Accurate Test Results
- Detection of Anti-HIV-1/2/O
- Detection of IgG and IgM Antibodies




ABON™ HIV Tri-Line Rapid Tests

A Preferred Choice in HIV Rapid Tests

Benefits

- Accurate Test Results
 - Sensitivity: 99.9%
- Specificity: 99.8%
- Detection of Anti-HIV-I(O) in Addition to Anti-HIV-1 and Anti-HIV-2
- Detection of IgG and IgM Antibodies
- Ideal for STAT Usage - Whole Blood Test Protocol & 10 Minutes Test Time
- Single Use Test's availability for STAT Usage

Test Procedure



Ordering information

Cat No	IHI-T402		IHI-T402H
Kit Size	40 Tests	25 Tests	10 Tests
Storage	2-30°C		
Sample Type	Whole Blood/Serum/Plasma		
Test Time	10 minutes		
Components	<ul style="list-style-type: none"> • Test devices • Droppers • Buffer • Package Insert 		<ul style="list-style-type: none"> • Test device • Disposable dropper • Buffer • Sterile lancet • Alcohol swab • Package insert



Inverness Medical Innovations Hong Kong Limited
 #15th Floor, Building B, Mordem International Plaza.
 159 Tianmu Shan Road, Hangzhou 310007, P.R.China
 Please visit our website for details: www.imihongkong.com
 E-mail: IMIHK.info@invmed.com

RAPID TEST ARGENTINA
 Rio Tercero
 Cordoba
 Argentina
www.rapidtest.com.ar
info@rapid-test.com.ar

Authorized Distributor

